



UT 2000 UNIVERSAL STRETCHER



Very light and versatile carrying device for the transport of injured persons in all kind of terrain. The UT 2000 can be divided in two halves for easy transport on the back of two rescuers. Each half can also be used as load carrying frame.

Tubular frame made from high-strength aluminum alloy. Special durable plastic lying area, shock and scratch-resistant, with good gliding abilities. Synthetic, rotting-free carrying-belts.

Wide and padded shoulderbelt with chest-fixation for optimum load distribution and high carrying comfort. Back netting for good air circulation. The carrying belts can be detached for helicopter transport. Four multi-functional load-bearing belts for the fixation of the injured or load.

Including **helicopter suspension belts**.

Load-carrying capacity: 2000 kg.



Art. No.	81301
Dimensions	Length: 1.810 mm / 2.000 mm Width: 440 mm High: 120 mm
Weight	Total: 8,0 kg
Accessories	One-Wheel System Art. 81320 Carrying Handles Art. 81322 Transport Bag Art. 81202



TYROMONT Alpin Technik GmbH
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria
T +43 (0) 5223 23 723-0 · office@tyromont.com
www.tyromont.com