



## UT 2000 ALPINE

### Universal Stretcher with integrated rescue bag



Rescue Bag open



Divided Universal Stretcher covering the Rescue Bag in one half of the stretcher

Further improvement of UT 2000 to better fit requirements of alpine rescue demands.

The integrated rescue bag allows an entire horizontal bedding of the injured person together with a vacuum mattress.

- The rescue bag is attached to the UT 2000 frame by a belt system on the bottom side of the bag.
- Quick fastener on one part of the rescue bag to divide the stretcher into two parts for transport.
- The wide area Velcro fastener allows an adjustable closing of the rescue bag.
- Strengthened head opening, geared up for the application of the optional Face Protection Shield (Art.No. 93281).
- The combined belt system fixates the injured person on the stretcher. Additionally the belt system can be used to transport both halves of the stretcher.
- Scope of delivery: UT 2000 Alpine with an integrated rescue bag, 4 fixation belts and helicopter suspension belts (Green carrying belt system of UT2000 not included)

#### Face Protection Shield (Art. No. 93281)

- Provides effective protection of the face and forehead area of the injured. Collapsible for easy and fast access. Transparent and cold-resistant foil material with Lexan battens.



Art. No.	<b>81200</b>	
Dimensions	Length:	2.000 mm
	Width:	440 mm
	High:	220 mm
Weight	Total (UT + Bag):	8,8 kg
	Total (Bag):	2,0 kg
Comments		



TYROMONT Alpin Technik GmbH

Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com

www.tyromont.com