

UT 2000 PACKTASCHE UT 2000 TRANSPORTATION BAG



XXL-Packtasche zur Aufbewahrung und Transport eines kompletten UT 2000 Sets.

Zur Aufbewahrung und Transport eines kompletten UT 2000 Alpin bzw. UT 2000 Sets. Verstaung der geteilten Rettungstrage, des Bergesacks (UT Alpin), zweier Führungsholme, zweier Radstützen und eines Rades. Das Rad wird dabei außen an der Packtasche befestigt.

Äußerst strapazierfähiges Material, inkl. ergonomisch geformte, gepolsterte und verstellbare Schultergurte mit Brustgurt für angenehmen Tragekomfort.

Geeignet für UT 2000 und UT 2000 Alpin mit Bergesack.

XXL transport bag for the transport and storage of a complete UT2000 set.

For the transport and storage of a complete UT2000 or UT2000 Alpin set. The transport bag offers place for the splitted stretcher, the rescue bag, two guiding handles, two wheelstruts and one wheel. The wheel gets fixed on the outside of the transport bag.

Made of really strong wearing material, with ergonomic formed, padded, and adjustable shoulder- and chestbelts for the best wearing comfort.

Suitable for UT2000 and UT2000 Alpin with rescue bag



P/N	81202	
Abmessung <i>Dimension</i>	Höhe / <i>Height</i> :	900 mm
	Breite / <i>Width</i> :	500 mm
	Tiefe/ <i>Depth</i> :	300 mm
Gewicht <i>Weight</i>	1,35 kg	
Anmerkungen <i>Comments</i>		