



SUSPENSION STANDARD for UT

Standard suspension for UT2000 and UT2000 Alpin.

Belts for a horizontal suspension consist of:

- Four black belts for head- and footarea with red loops.
- Two yellow belts for the lateral fixation with red loops.

The suspension length can be adjust due the red loops.

The suspension is included in the scope of delivery of the UT2000 (Art. 81301) and UT2000 Alpin (Art. 81200) with rescue bag.



Art. No.	81336	
Dimensions	Suspension length: 1.050 mm	
Weight	Total:	390 g
Comments		