



Vaccum Mattress 80 Mod. BAYERN with Cover and Carrying Hand-



Revolutionary new design: specially developed for the use with the **HRB TYROMONT Mod. BAYERN** (#93280).

Double-trapezoidal shape: provides extremely high side and length stability. The mattress remains very handy and light due to the cut-outs in the head and leg area.

Patented 48-Chamber-System: Holds the filling material always at its place and allows an extremely low profile of the vacuum mattress at maximum stability in vacuumized state. That's why the vacuum mattress can be instantly applied even in any terrain.

Stable inside-hull made from exceptionally resistant, fibre reinforced PU-material. Durable welds thanks to CAM weld-ing-procedure.

Low packing volume because of the little height of 50 mm of the mattress.

Rugged outside-cover provides extra protection for the mattress and is equipped with **eight carrying belts**. For clean-ing and disinfection, the cover could easily be removed.

Automatic, self-shutting valve does not have to be closed after vacuumizing the mattress. Even if the pump hose slips off, there is no loss of vacuum. The valve can be operated with any commercially available vacuum pump using an adapter.

Temperature resistance: -40°C to +70°C

Art. No.	81356	
Dimensions	Length:	2.000 mm
	Width:	800 mm
Weight	Total:	5,8 kg
Comments	Valve at head section.	



TYROMONT Alpin Technik GmbH
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria
T +43 (0) 5223 23 723-0 · office@tyromont.com
www.tyromont.com