



## VACUUM MATTRESS ULTRALIGHT with fixation belts and carrying handles



This new type of vacuum mattress offers several **advantages** over the previous vacuum mattresses:

**Light weight** only 4.7 kg.

**New multi-chamber system** - mattress maintains its shape. No more stroking necessary.

New **removable shoulder straps** and carrying straps.

Can be **rolled** up and easily transported with a shoulder strap.

The new **multi-chamber system (120 chambers)** - the layers of the outer cover are welded together point by point, creating chambers between the points that hold the inner beads in their correct position and are evenly distributed.

This enables **quick, safe and easy use** in all situations.

**Automatic, self-shutting valve** does not have to be closed after vacuumizing the mattress. Even if the pump hose slips off, there is no loss of vacuum. The valve can be operated with any commercially available vacuum pump using an adapter.

**Temperature resistance:** -40°C to +70°C



Art. No.	<b>81370</b>	
Dimensions	Length:	1.950 mm
	Width:	900 mm
Weight	Total:	4,7 kg
Comments	Valve at shoulder section.	



**TYROMONT Alpin Technik GmbH**  
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com  
[www.tyromont.com](http://www.tyromont.com)