



## ISOTHERM INLAY INJUREDBAG

The **Isotherm inlay** consists of two "blanket halves" that envelop the patient inside the mountain bag.

The inlays are hooked into the left and right sides of the mountain bag using Velcro. Can be **retrofitted** to existing bags.

This allows optimal use of the internal strapping while still providing a comfortable warm patient envelope.

This **bacteriostatic protective** cover keeps the patient warm thanks to the technical fiber HeatUP.

**Breathable, waterproof and fireproof** (M1) polyester cover coated with polyurethane, filling made of resistant fibers.

**Bacteriostatic:** viruses and bacteria cannot spread and multiply on the surface. Easy cleaning with quick disinfectant wipes and ready to use again in a few minutes.



Art. No.	<b>92105</b>	
Dimensions	Length:	1.750 mm
	Width:	700 mm
Weight	Total:	1,6 kg
Comments	<b>retrofitted</b> to existing bags	