

ISOTHERM-MANTEL ISOTHERM-COAT



Leichtes und kompaktes System zur einfachen Bergung von unverletzten, jedoch erschöpften oder unterkühlten Personen, wie z. B. Trailrunnern im alpinen Gelände. Unterstützt den Wärmeerhalt und reduziert eine schnelle Hypothermie. Das modulare Konzept ermöglicht eine Anpassung an unterschiedliche Wetter- und Einsatzbedingungen. Durch das geringe Packmaß ist das System leicht transportierbar und platzsparend verstaubar.

Der verwendete Stoff ist bakteriostatisch, reduziert das Wachstum von Bakterien und Viren und lässt sich **einfach mit handelsüblichen Desinfektionstüchern reinigen**.

Lightweight and compact system designed for the easy evacuation of uninjured but exhausted or hypothermic persons, such as trail runners in alpine terrain. It supports heat retention and helps reduce rapid hypothermia.

The modular concept allows adaptation to different weather and operational conditions. Its small pack size makes the system easy to transport and space-saving to store.

The material is bacteriostatic, reducing the growth of bacteria and viruses, and can be easily **cleaned using standard disinfectant wipes**.

P/N	82110	
Abmessung Dimension	Länge / <i>Length</i> :	1000 mm
	Breite / <i>Width</i> :	500 mm
Gewicht Weight	0,6 kg	
Anmerkungen Comments		