



HELICOPTER RESCUE BAG TYROMONT Mod. TACTICAL with Internal Patient Restraint System



HEMS (=Helicopter Emergency Medical System) carrying bag for the horizontal transport of injured persons.

Protective bag to fully enclose the patient, manufactured from high-strength, water- and wind resistant perlon PES fibre.

Bottom with **sewn diagonal-belt weight bearing system** made from rotting-free, high strength 45mm PES belt provides even weight distribution for perfectly horizontal bedding of the injured.

Underside equipped with **spike-stripes**. Integrated **carrying loops** for lifting the patient and attachment of the suspension, **colour coded** for the attachment of an optional detachable suspension.

Stow-away side-pouches over the entire length for the complete storage of the suspension.

Longitudinal large-area Velcro closing for the perfect adoption to the size of the injured. Four additional external **fixation belts**.

Integrated vacuum suction system for the operation of the vacuum mattress through an external suction point.

Equipped with an **internal patient restraint system (=IPRS)** for the securing of the injured in longitudinal and lateral direction, additionally equipped with leg loops for fixation of the injured in caudally direction (certified according EN 1497:2007).

Four **retaining tabs connected with high-load Cobra® buckles** for the fixation of the bag with the secured injured in the cabin at its hard points.

Art. No.	922X - 92214	
Dimensions	Length:	2.050 mm
	Width:	560 mm
Weight	Total:	3,7 kg
Reccomended Accessories	Transportation Bag TACTICAL	(Art. 98014)
	Face Protection Shield	(Art. 93281)
	Vacuum Mattress 80 TAC.	(Art. 81357)
	Vacuum Pump EVO	(Art. 93134)
	Rotation Brake System TAC.	(Art. 93354)
Certification	EC type examination no. E 7111 (NB 1246) (according EASA CM-CS-005, Simple PCDS).	



TYROMONT Alpin Technik GmbH
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria
T +43 (0) 5223 23 723-0 · office@tyromont.com
www.tyromont.com