



SPECIAL RESCUE BAG - TRB



Rescue bag with comfortable dampening laying area designed for the usage without other mats.

Unsurpassed shelter against cold, snow and wetness and effective protection against ground impacts for the injured.

Thick, comfortable dampening mat plus temperature insulation layers at the bottom and both sides.

Completely manufactured from **high-strength, water-resistant nylon-fabric**. This high-grade material ensures longevity, maximum wear-resistance, pleasant handling and does not become stiff even under very cold temperatures. Additional reinforcements on the bottom guarantee years of use.

Strong, freeze-safe, wide area Velcro closing, allows tight closing and adjustment to the injured.

Padded hood to cover the head.

The **foot compartment** can be closed separately and provides enough volume to accommodate even large ski boots and/or splints.

Eight carrying loops for comfortable carrying of the injured.

Suitable for all TYROMONT rescue sleds or other models.



Art. No.	93008	
Dimensions	Length:	2.120 mm
	Width:	560 mm
Weight	Total:	3,9 kg
Comments	With comfortable dampening mat.	

