



## AKJA 2200 PLUS



Designed for **2-rescuer** operation in **all kind of snow** - covered terrain, groomed or ungroomed.

Especially good in **steep and/or icy slopes** and moguls.

The **classic Akja - Rescue Sled** in one-piece version. Timeproofed design and workmanship, thousands of them worldwide running.

+ 20 cm in length - more room for even the tallest injured  
+ 20 cm in height of the poles for more comfort of the rescuer

**Long lasting durability:** Manufactured completely from high-strength, corrosion-resistant aluminum alloy for many years of reliable, maintenance-free action.

**Stable running characteristics:** Reinforced steel runners with especially hardened fins provide stable running even in very icy conditions.

**Two pairs of poles** (for two rescuers), can be detached easily by bayonett-lockings. The grips of the handles are high enough for tall rescuers. With integrated stabilization profiles for increased bending strength.

**Bayonett-lockings** are wear-free, 100% icing-safe and easily to handle.

**Carrying grips** on both ends of the sledge for comfortable lifting and carrying of the sled.

**Four fixation belts with aluminum buckles** for safe fixation of the injured; resistant against breaking and icing, can be operated with skiing-gloves.

Art. No.	<b>93031</b>	
Dimensions	Length:	2.200 mm
	Width:	560 mm
	Length of poles:	1.650 mm
Weight	Total: (incl. Poles)	27,0 kg
Comments	On request length of 2.000mm available (Art. 93059)	



**TYROMONT Alpin Technik GmbH**  
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com  
[www.tyromont.com](http://www.tyromont.com)