



## AKJA 2000 DIVISIBLE



**Now** available in the **standard length of 220cm**.

- **Robust dividing mechanism** guarantee stable driving characteristics even in rough terrain
- Divisible in two halves for **easy transport**
- The length of 220cm (110cm divided) represents a good compromise between lying comfort for the injured and carrying comfort when transporting the rescuer.

**Akja for two-man operation** for ski areas in all difficulties, especially for unprepared and / or impassable terrain.

**Years of durability:** The slide body is made entirely of high-strength aluminum alloy. All steel parts are either galvanized or powder-coated.

**Very robust dividing mechanism.** Stable driving characteristics: steel runners with crampons on the underside of the Akja and additional lateral runners for safety on all types of snow.

Four **removable drawbars** with large round handles for a safe and comfortable hold. Particularly resistant to bending thanks to the integrated reinforcement profile. Also available with divisible drawbars on request.

**Drawbar bracket with bayonet quick-release** fasteners, wear-free and robust. Easy and safe to use, 100% ice-proof.

**Carrying handle** at both ends for comfortable lifting and carrying of the Akja with the injured person.

**Four restraint belts with aluminum clamping buckles** for a secure hold of the injured person. Can be operated with

Art. No.	<b>93036</b>	
Dimensions	Length:	2200 mm
	Width:	560 mm
Weight	Total incl. drawbar:	28 kg
Comments		