



AKJA 2150 COMBI

New development for most versatile operation.

Based on the design of the well-proved standard Akja, this sled is tailored for both **single-rescuer operation** in easy to medium terrain and **two-rescuer operation** in difficult terrain.

New, **flat stern, longer runners** on the bottom and **special fins** to keep it stable in single-rescuer operation.

Additionally, the new, even stern allows a **perfectly horizontal bedding of the injured** at a very low centre of gravity (essential for riding-stability).

Long lasting durability: Manufactured completely from high-strength, corrosion-resistant aluminum alloy for many years of reliable, maintenance-free action.

Two pairs of poles (for two rescuers), can be detached easily by bayonett-lockings. The grips of the handles are high enough for tall rescuers. With integrated stabilization profiles for increased bending strength.

The **second pair of handles** can be **attached at the side** of the sled during single-rescuer operation.

Bayonett-lockings are wear-free, 100% icing-safe and easily to handle.

Carrying grips on both ends of the sledge for comfortable lifting and carrying of the sled.

Four fixation belts with aluminum buckles for safe fixation of the injured; resistant against breaking and icing, can be operated with skiing-gloves.



Art. No.	93047	
Dimensions	Length:	2.150 mm
	Width:	560 mm
	Length of poles:	1.650 mm
Weight	Total incl. poles:	28,0 kg
Comments	Suitable for single-rescuer operation.	

