



THERMO RESCUE BAG—TRB



Newly developed rescue bag for the ideal transport of injured in alpine terrain.

High-quality rescue bag for use with Akja, Canadier, Mountain Stretcher or any other type of rescue device.

Ergonomic design, **high-strength, water-resistant nylon-fabric**, with strong wide area Velcro fastener for adjustable closing of the bag.

Two additional zippers at the upper part of the bag allow a fast access to the injured person (e.g. for an infusion).

Foot-compartment with two zippers can be opened separately and provides enough space for large ski boots and/or splints.

Thermo insulation layers at both sides and the bottom that can be removed for cleaning purposes.

Six tight carrying loops for comfortable carrying of the injured.

Very **compact foldable**, little weight of 2,0 kg.

Without additional dampening mat at the bottom; therefore, optimum application in **combination with a vacuum mattress**.

Suitable for all TYROMONT rescue sleds, mountain stretchers and other models.



Art. No.	93070	
Dimensions	Length:	2.050 mm
	Width:	550 mm
Weight	Total:	2,0 kg
Comments	Ideal combination with a vacuum mattress.	



TYROMONT Alpin Technik GmbH
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria
T +43 (0) 5223 23 723-0 · office@tyromont.com
www.tyromont.com