



PROTECTION BAG FOR STRETCHER



Provides **optimum stabilization and entire weather protection** for the injured while **transported on the stretcher** of the TYROMONT rescue trailer.

The **purpose-built bag** for the stretcher is solidly mounted on the stretcher and provides, together with its large Velcro fastener, stable hold for the injured, even in rough terrain. The bag is wide enough to accommodate additional splints and blankets.

Thanks to its special nylon-fabric, the bag provides **effective wind and water protection** and could be easily cleaned.



Art. No.	93076	
Dimensions	Length:	1.950 mm
	Width:	500 mm
Weight	Total:	1,8 kg
Accessories	Back Cover for Rescue Trailer	Art. No. 93164
	Rescue Trailer Tyromont	Art. No. 93152
	Stretcher Rescue Trailer	Art. No. 93162
	Vacuum PAD 60 - STD	Art. No. 81354
	Vacuum Mattress 80 - STD	Art. No. 81350
	Vacuum Pump „EVO“ with foot bow	Art. No. 93134