



PATIENT BAG FOR BUCHER/ FERNO STRETCHER

Developed for various helicopter carrying systems.

TYROMONT has developed a special patient bag in order to **be able to transport the patient safely** from external influences even during helicopter transport on various different stretchers (e.g. Bucher, Stryker or remote stretcher).

This system offers both the crew and the patient a multitude of different options.

- # **Handles** for the patient
- # **Spit pocket** in the head area
- # **Fully openable/removable coverlet** - provides the clinician with easy patient management and patient thermal management.

Vacuum mattress/immobilization products can be used/ integrated as usual. **Attachment loops** to attach the patient bag to the stretcher.

Easy removal of the bag for **cleaning**. Cleaning with a disinfectant wipe or washing machine possible. Pictures show different patient bags.

The **patient bags** can be developed and modified **individually** depending on the **stretcher** or **customer requirements**.



Schlitz zur
Befestigung des
Aufsatzfisches



Patientensack geöffnet

| | | |
|------------|--------------|----------|
| Art. No. | 93082 | |
| Dimensions | Length: | 1.950 mm |
| | Width: | 500 mm |
| Weight | Total: | 4,4 kg |
| Comments | | |