



## INJURED BAG TYRAL



The TYRAL injury protection sack **completely encloses** the patient (also with a vacuum mattress) and protects against impairment and hypothermia.

An **ice-proof** Velcro fastener enables injured persons to be individually sized.

The **head opening** of the rescue bag is equipped with a reinforced edge and contains the holder for the optional face protection hood (item no. 93281).

Equipped with an **A5 document pocket** with a large viewing window.

**Eight carrying straps** that can be stowed in the side pockets allow the injured person to be lifted or repositioned.

With **integrated restraint belt** (= IRG) to secure the patient lengthways and crossways. Additional leg loops to fix the injured person in the caudal direction prevent the injured person from slipping in the direction of the foot. Four external COBRA® buckles for securing the sack to the optional fastening loops of the TYRAL mountain stretcher (item no. 70550).



|            |                                                               |
|------------|---------------------------------------------------------------|
| Art. No.   | <b>93084</b>                                                  |
| Dimensions | Length: 2.050 mm<br>Width: 450 mm<br>Coverage: 1.400—1.700 mm |
| Weight     | Total: 3 kg                                                   |
| Comments   |                                                               |

