



## SALVAGETRANSPORTBAG 7 - TERRA



The **salvage transport bag 7-Terra** has been designed for the transport of fatally injured persons in the alpine terrain.

It is characterized by its light weight and easy handling. The **12 carrying loops** make it possible to carry the casualty safely in pathless terrain.

The large zipper on three sides allows easy opening and closing of the transport bag at any time.

On the front side there is an **A5 document pocket** with a large viewing window for storing important information.

The **welded seams** inside the bag prevent water and other liquids from escaping unintentionally.

The use of tarpaulin fabric allows **easy and quick cleaning** with water and a light soap suds without long drying time.



Art. No.	<b>93090</b>	
Dimensions	Length:	2.020 mm
	Width:	500 mm
Weight	Total:	3,0 kg
Comments	Supplied with transport bag	



**TYROMONT Alpin Technik GmbH**  
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com  
[www.tyromont.com](http://www.tyromont.com)