



## VACUUM LEGSPLINT - STD

Vacuum product for the optimum splinting at all leg injuries.

**For the optimum splinting at all leg injuries** up to the femur. Due to belt-system, the legsplint can be adjusted exactly to the leg before vacuumizing. Therefore, the rescuer has both hands free during vacuumizing.

**Very stable, fibre-reinforced PU-cover**, multichamber system for even and non-skidding distribution of the filling material.

**Automatic, self-shutting valve.**



Art. No.	<b>93107</b>	
Dimensions	Length:	1.100 mm
	Width:	400 - 700 mm
Weight	Total:	1,6 kg
Comments		