



## VACUUM ARMSPLINT STD



Vacuum product for the optimum splinting at all arm injuries.

**For the optimum splinting at all arm injuries.** Due to belt-system, the armsplint can be adjusted exactly to the arm before vacuumizing. Therefore, the rescuer has both hands free during vacuumizing.

**Very stable, fibre-reinforced PU-cover,** multichamber system for even and non-skidding distribution of the filling material.

**Automatic, self-shutting valve.**



Art. No.	<b>93108</b>	
Dimensions	Length:	670 mm
	Width:	290 - 480 mm
Weight	Total:	0,7 kg
Comments		



**TYROMONT Alpin Technik GmbH**  
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com  
[www.tyromont.com](http://www.tyromont.com)