

VAKUUMPUMPE EVO mit Fussbügel VACUUM PUMP EVO with Foot Bow



Hocheffiziente Aluminium-Vakuumpumpe zum raschen Absaugen von Vakuummatratzen.

Stabile Vakuumpumpe mit Fußbügel aus Aluminium zum Absaugen der TYROMONT Vakuummatratzen und Schienen. Integriertes Rückschlagventil sorgt für ein halten des Vakuums, wenn die Pumpe losgelassen wird.

Mit kältebeständigem Silikonschlauch und neu entwickeltem Metalladapter zum raschen und hocheffizienten Absaugen von Vakuummatratzen und –schienen.

Abgewinkelte Schlauchführung für optimalen Schutz des Schlauchs vor Beschädigungen durch Ski- oder Bergschuhe.

Durable pump with foot bow for use with all TYROMONT vacuum mattresses and splints.

Robust vacuum pump with aluminium foot stirrup for suctioning TYROMONT vacuum mattresses and splints.

Integrated **non-return valve** ensures that the vacuum is maintained when the pump is released.

With **cold-resistant silicone tube** and newly developed metal adapter for quick and highly efficient suctioning of vacuum mattresses and splints.

Angled tube guide for optimal protection of the tube from damage by ski or mountain boots.



P/N	93134	
Abmessung Dimension	Länge Pumpe/ Length pump:	700 mm
	Länge Schlauch/ Length Hose:	525 mm
Gewicht Weight	1 kg	
Anmerkungen Comments		