



## RESCUE TRAILER TYROMONT for snow-mobiles and quads



Injured stretcher fixed on trailer.



Guiding runners with fins.



Ground anchors for fixation of stretcher on a slope.



**High-grade, completely equipped rescue trailer for snow-mobiles and quads.**

Optimum protection of the injured against snow-spray, wind and cold.

**Maximum safety** by two rollover bars and fixation belts with aluminum buckles.

**Integrated, extractable stretcher** with adjustable reclining area with thermo-insulation. Stretcher equipped with diagonal fixation belts for the upper part of the body and transverse belts for the hip and legs.

Additionally equipped with two **ground anchors** for fixation of the stretcher on snow or frosty underground.

Guiding profiles and **quick-locks** for easy loading of the injured on the stretcher.

Additional **guiding profiles** for optional loading of an Akja or Canadier rescue sledge (please specify your demand during ordering).

**Cover** manufactured from water- and wind-resistant PES fibre with large, cold-resistant clear-areas.

Cover could be opened by a **large, framed zipper** to easily handle the injured person.

**Stable running properties** through wide track, guiding runners and additional four fins.

**Dampened suspension** for a smooth ride of the injured.

For **summer operation** convertible to **wheels**.

Art. No.	<b>93152</b>
Dimensions	Total: 2.250 x 790 mm Length trailer: 2.850 mm Width incl. runners: 1.080 mm Height incl. Cover: 1.450 mm
Weight	Total incl. rescue set: 105 kg
Comments	Use with Canadier and injured stretcher After-market equipment without trailer



**TYROMONT Alpin Technik GmbH**  
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com  
[www.tyromont.com](http://www.tyromont.com)