



SUSPENSION PRO TYRAL



Suspension Pro TYRAL mountain stretcher.

Sling leg area 120 cm; Orange 21mm
Loop head area 100 cm; Yellow 21mm

Attachment of the stretcher using a Twistlock carabiner
Attachment to the rope on 2 rigger plates and carabiners.

The suspension can be converted from a horizontal to a vertical suspension without tools.
Simply unthread the two yellow ribbon loops from the rigging plate and use them for vertical lowering.

Art. No.	93257	
Dimensions	Length Head:	100 mm
	Length Feetarea:	120 mm
Weight	Total:	0,8 kg
Comments	Delivered with TYRbag	

