



## COBRA CONNECTION SET FOR CARRYIN STRAPS TYRAL



This kit includes four COBRA mounts designed for retrofitting four carrying straps to the GT TYRAL mountain stretcher. The mounts can be quickly and easily hooked into the designated slot between the frame tube and the oval handle—no tools required. This ensures fast, secure, and flexible transport of the stretcher, making it ideal for use in rough terrain or rescue operations.

### Scope of delivery:

4 COBRA mounts for use with 4 carrying straps.



Art. No.	93259	
Weight	Total:	0,22 kg per pair
Comments		



**TYROMONT Alpin Technik GmbH**  
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria  
T +43 (0) 5223 23 723-0 · office@tyromont.com  
[www.tyromont.com](http://www.tyromont.com)