



EVACUATION BOARD TYRAH N-E-xt



Further development of the TYRAH with spinal stabilization.

By using the special Tyrah rescue triangle, the NExt guarantees the immobilization and support of the spine during rope maneuvers or rescue with rope technique.

The special rescue system with rope or winch maneuver ensures the best position of the patient during handling and avoids problems related to compartment syndrome.



Art. No.	93268	
Dimensions	Length:	900 mm
	Width:	320 mm
Weight	Total:	2,5 kg
Accessories	Rescue Seat TYRAH AR (93222) TYROLL CS (93265)	

