



FACE PROTECTION SHIELD



Foldable and removable protective shield for helicopter rescue bags, made from extremely durable, transparent special film with integrated PE3 strips.

The hood protects the **patient's facial area** from small stones, dust, ice, snow, and helicopter downwash. Its open design allows for continuous contact with the patient.

A **light-proof version** of the protective hood is also available and is used during night flights, where the patient may be blinded by the helicopter's searchlights.

The face protection hood is compatible with the following products:

- HBS TYROMONT Model BAYERN
- HBS TYROMONT Model CHRISTOPHORUS II
- HBS TYROMONT Model CHRISTOPHORUS EVO
- Rescue bag for UT2000 Alpine
- Injury Protection Bag GT TYRAL



Art. No.	93281 / 93349	
Variants	93281 – Transparent Shield 93349 – Lightproof Shield	
Weight	Total:	0,3 kg
Comments		