



RESCUE SLING TYRloop

The **Rescue Sling TYRloop** is used for rescuing people in water and land emergency situations.

Due to the **one-size closure system** (Velcro/hook closure of the belt), the adjustment to the size of the person to be rescued is automatic or self-adjusting, and the patient's own weight ensures that it is securely wrapped around the patient.

The **integrated flotation device** with three external hand-holds allows rescue from water emergency situations. **This does not serve as a buoyancy aid** for the patient, it only allows the rescue loop to float up.

Take care to avoid suspension trauma of the person to be rescued.



Art. No.	93359	
Dimensions	Length total:	2.200 mm
	Length loop:	900 mm
Weight	Total:	0,7 kg
Comments	EN 1498 class A EASA CM-CS-005 Simple PCDS	