



HOR-VER_SYSTEM II WITH GRILLON 5m



System for horizontal and vertical lowering of the rescue bag.

The HOR-VER-Sys II suspension system is used to rescue people from **confined spaces**. The system can be positioned from the horizontal to the vertical position by one person without much effort due to the **integrated pulley**.

The patient can be safely brought back to a horizontal position by means of a grillon on the pulley.

This system may only be used with rescue bags in combination with an **integrated restraint belt** to secure the patient in the longitudinal and transverse direction and four external Cobra buckles.

The HOR-VER-Sys II is supplied with a matching bag.



Art. Nr.	93534	
Gewicht	Gesamt:	2,3 kg
Anmerkungen		