



## 2F-ADJUSTABLE SUSPENSION VERTICAL

Adjustable vertical suspension for TYRAL mountain stretcher.

Using the **rigger plate** as an **anchor point**, the **two adjustable slings** become a simple tool for vertical abseiling situations.

The adjustability offers the **ideal length** for all types of terrain to ensure a simple and **smooth abseiling process**.

**Length adjustable from 30cm to 110cm.**

Attached to the stretcher using a **twistlock carabiner**.



Art. No.	<b>93567</b>	
Dimensions	Length:	30 cm - 110 cm
Weight	Total:	1 kg
Comments		