



TRANSPORTBAG TYRAL AIR



The **TYRAL AIR transport bag** is intended for use in bad weather.

The TYRAL mountain stretcher can be flown up to the cloud cover and dropped off there with the mountain rescuers.

All-round straps and carrying handles.

Two carabiners for attaching to the helicopter rope.

Helicopter suspension can be **stowed** away.

Stable bag with a payload of **up to 120 kg**.



Art. No.	93735	
Dimensions	Length:	1.200 mm
	Depth:	350 mm
	Width:	800 mm
Max. load	120 kg	
Comments		