



WINTERMODUL TYRAL

The **winter module** for the TYRAL mountain stretcher consists of a **monoski and extended bars**.

The monoski is mounted over the front edge of the stretcher and, like the two sliding ski modules, can be **easily attached without tools**.

The two sliding ski modules result in **driving characteristics** that no mountain skier has ever achieved in the snow. (please order sliding ski module with nr.:93262)

The two extended bars can be **easily divided and easily transported** in three parts (in or on the backpack).

These are mounted on the existing bars without tools and offer stable properties when driving the mountain stretcher in the outback.

This makes the **TYRAL** mountain stretcher suitable for **use in all 4 seasons** and one of the lightest and most universally applicable alpine rescue devices.



Art. No.	93741	
Material	Aluminium	
Weight	Total:	app. 5 kg
Comments		