



DOUBLE WHEEL SYSTEM GT TYRAL



Double Wheel for the TYRAL Mountain Stretcher

The double wheel for the TYRAL mountain stretcher is a robust and reliable solution designed specifically for versatile applications. The wheel support is attached to the stretcher, as before, using two quick-release pins, enabling quick and secure installation.

Features:

- **Sturdy Construction:** Two large tires provide excellent stability for the stretcher, making it ideal for various terrains, such as mountain trails, forest paths, as well as narrow alleys, campsites, and festivals.
- **Tire Dimensions:** Each tire has a diameter of 400 mm and a width of 100 mm, ensuring optimal handling and high load capacity.
- **Excellent Rolling Properties:** Thanks to the built-in ball bearings in the tires, the double wheel offers outstanding rolling properties, making it easier to maneuver even on challenging surfaces.

Applications:

The double wheel is versatile and well-suited for:

- Mountain and forest trails
- Narrow alleys
- Campsites
- Festivals
- Well-maintained hiking paths

Art. No.	93743	
Material	Wheel support made of high-strength aluminium	
Weight	Total:	9,8 kg
Comments		

