



E-WHEEL GEBIRGSTRAGE TYRAL



Electric motor assisted wheel for easy transport of equipment to the scene of the accident. The electric motor also has a braking effect when transporting the patient downwards.

Wheel with powerful electric motor for the TYRAL mountain stretcher.

Easily transport equipment to the scene of the accident.

The electric motor **dims the speed** using a remote control.

The **remote control** can be operated manually by a third person and thus enables optimal transport. When traveling downwards, the motor can also serve as a brake.

The two batteries are long-lasting and provide power for around 2 hours each at normal temperatures. The **battery module is mounted on the wheel support**.

The **battery pack can be replaced** without tools while driving.

Scope of delivery: Wheel with electric motor, 2x battery pack and wheel support attachment, remote control.



Art. No.	93747
Weight	Total (without wheel support): 8 kg
Comments	



TYROMONT Alpin Technik GmbH
Bert-Köllensperger-Str. 6 · 6065 Thaur/Austria
T +43 (0) 5223 23 723-0 · office@tyromont.com
www.tyromont.com