



HEAD PROTECTION FELX



The lightweight rockfall protection for the head area has been specifically developed for use with the **TYRAL mountain rescue stretcher** and provides additional protection for the patient during alpine rescue operations. Especially during vertical evacuations or lowering operations, the head protection reliably shields the injured person from falling rocks, ice, or other objects.

The protection can be easily attached to the stretcher frame using Velcro straps, allowing for quick installation or removal. The integrated shock-absorbing material reduces impact forces and is also used in high-quality back protectors in the ski industry, ensuring a high level of protection while maintaining a very low weight.

For transport, the rockfall protection can be rolled up compactly and easily stored in a backpack, ensuring it is always readily available without taking up unnecessary space. The head protection is fully compatible with both the adjustable and the fixed vertical suspension systems from TYROMONT and, in addition to the **TYRAL mountain stretcher**, is also suitable for use with the **STEELIGHT mountain rescue stretcher**.



Art. No.	93767	
Dimensions	Length:	500 mm
	Width:	630 mm
Weight	Total:	0,6 kg
Comments		

