

VAKUUMPUMPEN ADAPTER ADAPTER FOR VACUUM PUMP EVO

Metall-Ersatzadapter für Vakuumpumpen zum hocheffizienten Absaugen der TYROMONT Vakuummatratzen.

Dieser innovative, **L-förmige Metalladapter** wurde speziell für das **effiziente Absaugen** von TYROMONT Vakuummatratzen und Vakuumschienen entwickelt. Durch die präzise Verarbeitung und die hochwertige Materialwahl garantiert er eine zuverlässige, dichte Verbindung während des gesamten Absaugvorgangs.

Ausgestattet mit integrierten Abdichtungs- und Fixierungsringen sorgt der Adapter für eine sichere und luftdichte Verbindung – **auch unter anspruchsvollen Einsatzbedingungen**. Die Passgenauigkeit für das selbstschließende TYROMONT Vakuum-Ventil ermöglicht eine schnelle, einfache und verlustfreie Handhabung.

Geeignet für Schlauchdurchmesser von 6 bis 10 mm – ideal für den professionellen Rettungs- und Sanitätseinsatz.



Metal- reserve adapter for vacuum pumps for a high efficient vacuumizing of TYROMONT vacuum mattresses and splints.

Newly developed, **L formed metal adapter** with lining- and fixation rings for a **high efficient vacuumizing** of TYROMONT vacuum mattresses and splints.

Accurately fitting to the self-shutting TYROMONT vacuum valve.

P/N	93137
Gewicht <i>Weight</i>	24 kg
Anmerkungen <i>Comments</i>	